




SNAP-OFF

SCREW SYSTEM





SNAP-OFF

SCREW SYSTEM

Quick - Low-Profile - Superior Purchase



STEP 1 - Determine Screw Length


Measure the desired screw length with a ruler to the plantar side of the osteotomy.

STEP 2 - Advance Screw

Advance the screw while securing the osteotomy in its correct orientation.

STEP 3 - Snap Off Screw Head

Ensure the Pin Driver is fully seated to the base of the Snap-off Screw head. Bend the screw shaft with the Pin Driver until the wire breaks away from the screw head.



Indicated for use in bone reconstruction, osteotomies, arthrodesis, joint fusion, fracture repair, & fracture fixation.

Product Number	Diameter	Length
M10-20411	2.0mm	11mm
M10-20412	2.0mm	12mm
M10-20413	2.0mm	13mm
M10-20414	2.0mm	14mm
M10-20415	2.0mm	15mm
M10-20416	2.0mm	16mm
M10-20417	2.0mm	17mm
M10-20418	2.0mm	18mm
M10-20419	2.0mm	19mm
M10-20420	2.0mm	20mm

Fuse Medical, Inc.

4343 Sigma Road, Suite 500
Farmers Branch, TX 75244
469.862.3030

info@fusemedical.com
www.fusemedical.com

SNOF-B-001 REV C Effective date: 09/23/2022 DC-0427

Compression Screw System