



TyWEDGE™ SYSTEM

EVANS AND COTTON PROCEDURES

TYPEEK™ TITANIUM PLASMA SPRAY





PEEK RADIOLUCENT BODY

BIOACTIVE SURFACES


OPTIMIZED TOOTH PROFILE

BONE GRAFT CONTAINMENT

COMPLETE INSTRUMENTATION





MADE IN AMERICA



*Titanium plasma sprayed coating on PEEK has been shown to encourage osseointegration [Walsh et. al., 2015]



Implants - Evans



Height	Length x Width		
	18 x 18mm	20 x 20mm	22 x 22mm
8	X06-E01808	X06-E02008	X06-E02208
10	X06-E01810	X06-E02010	X06-E02210
12	X06-E01812	X06-E02012	X06-E02212

Height	Approximate Autograft Volume		
	18 x 18mm	20 x 20mm	22 x 22mm
8	0.51 cc	0.74 cc	1.00 cc
10	0.59 cc	0.85 cc	1.15 cc
12	0.67 cc	0.96 cc	1.31 cc

Implants - Cotton



Height	Length x Width	
	16 x 14mm	20 x 14mm
4.5	X06-C01604	X06-C02004
5.5	X06-C01605	X06-C02005
6.5	X06-C01606	X06-C02006

Height	Approximate Autograft Volume	
	16 x 14mm	20 x 14mm
4.5	0.17 cc	0.25 cc
5.5	0.20 cc	0.27 cc
6.5	0.22 cc	0.31 cc

Corporate Headquarters

Fuse Medical Inc.

1565 N Central Expressway Suite 220
Richardson, TX 75080

P. 469.862.3030 | F. 469.862.3035 | FUSEMEDICAL.COM | OTC: FZMD