



STERIZO™

TOTAL KNEE SYSTEM

STABILITY
STRENGTH
SIMPLICITY



In a HIGH FLEX DESIGN



Stability

Choose your constraint style – Cruciate Retained, Ultra-Congruent/ Anterior Stabilized, or Posterior Stabilized, the Sterizo Total Knee System provides stability for your patients throughout the entire range of motion.

Strength

Sterizo literally means strong. The Sterizo Total Knee System is made in America, bringing you strength in confidence. Strength in quality. Strength in longevity. Strength in patient outcomes and satisfaction.

Simplicity

The Sterizo Total Knee System was designed with efficiencies in mind – both in elegant implant designs to accommodate all patient anatomies, and streamlined, modular instrumentation – all to assist you with simplifying the surgical procedure and improving Operating Room utilization.


High Flex Design

The Sterizo Total Knee System was engineered to deliver high flexion across all constraint platforms. 139° of flexion allows your patients to achieve strong, reproducible results, all while maintaining anatomic high-flex mobility.

STERIZO™


COMPONENT DIMENSIONS

Femoral Component - CoCrMo Alloy



	ML	AP
	(MM)	(MM)
Sz.1	60	56
Sz.2	64	60
Sz.3	68	64
Sz.4	72	68
Sz.5	76	72
Sz.6	80	76

Tibial Tray - Titanium Alloy




	ML	AP
	(MM)	(MM)
Sz.1	66	44.5
Sz.2	69	47
Sz.3	72	49.5
Sz.4	76	52.5
Sz.5	80	55.5
Sz.6	84	58.5

Tibial Insert - UHMPE - Highly Cross Linked



THICKNESS
(MM)
9
11
13
15
18

Onset Patella - UHMPE



THICKNESS	DIAMETER
(MM)	(MM)
8	29
8	32
9	35
9	38
10	41

ABOUT FUSE MEDICAL :

Fuse is an emerging manufacturer and distributor of innovative medical devices for the orthopedic and spine marketplace. We provide a comprehensive portfolio of products in the orthopedic total joints, sports medicine, trauma, foot and ankle, as well as, degenerative and deformity spine, orthobiologics and regenerative medicine products.



1565 N. Central Expressway
Suite 220
Richardson, TX 75080
P: 469-862-3030
F: 469-862-3035
info@fusemedical.com
www.FUSEMEDICAL.com
OTC: FZMD