



Compression Screw Set Cannulated - Headed



- **Titanium**
- **Self drilling**
- **Self tapping**
- **PickNStick**

MADE IN AMERICA



Compression Screw Set Cannulated - Headed

Length (mm)	2.5mm Headed Cannulated Screw	3.0mm Headed Cannulated Screw	4.0mm Headed Cannulated Screw		Length (mm)	4.5mm Headed Cannulated Screw		6.5mm Headed Cannulated Screw		
	Long Thread	Long Thread	Short Thread	Long Thread		Short Thread	Long Thread	Short Thread	Long Thread	Full Thread
10	M9-25310	M9-30310			30	M9-45230	M9-45330	M9-65230	M9-65330	M9-65430
11	M9-25311	M9-30311			32	M9-45232	M9-45332	M9-65232	M9-65332	M9-65432
12	M9-25312	M9-30312			34	M9-45234	M9-45334	M9-65234	M9-65334	M9-65434
13	M9-25313	M9-30313			36	M9-45236	M9-45336	M9-65236	M9-65336	M9-65436
14	M9-25314	M9-30314	M9-40214	M9-40314	38	M9-45238	M9-45338	M9-65238	M9-65338	M9-65438
15	M9-25315	M9-30315	M9-40215	M9-40315	40	M9-45240	M9-45340	M9-65240	M9-65340	M9-65440
16	M9-25316	M9-30316	M9-40216	M9-40316	42	M9-45242	M9-45342	M9-65242	M9-65342	M9-65442
17	M9-25317	M9-30317	M9-40217	M9-40317	44	M9-45244	M9-45344	M9-65244	M9-65344	M9-65444
18	M9-25318	M9-30318	M9-40218	M9-40318	46	M9-45246	M9-45346	M9-65246	M9-65346	M9-65446
19	M9-25319	M9-30319	M9-40219	M9-40319	48	M9-45248	M9-45348	M9-65248	M9-65348	M9-65448
20	M9-25320	M9-30320	M9-40220	M9-40320	50	M9-45250	M9-45350	M9-65250	M9-65350	M9-65450
21	M9-25321	M9-30321	M9-40221	M9-40321	52	M9-45252	M9-45352	M9-65252	M9-65352	M9-65452
22	M9-25322	M9-30322	M9-40222	M9-40322	54	M9-45254	M9-45354	M9-65254	M9-65354	M9-65454
23	M9-25323	M9-30323	M9-40223	M9-40323	56	M9-45256	M9-45356	M9-65256	M9-65356	M9-65456
24	M9-25324	M9-30324	M9-40224	M9-40324	58	M9-45258	M9-45358	M9-65258	M9-65358	M9-65458
25	M9-25325	M9-30325	M9-40225	M9-40325	60	M9-45260	M9-45360	M9-65260	M9-65360	M9-65460
26	M9-25326	M9-30326	M9-40226	M9-40326	65	M9-45265	M9-45365	M9-65265	M9-65365	M9-65465
27	M9-25327	M9-30327	M9-40227	M9-40327	70	M9-45270	M9-45370	M9-65270	M9-65370	M9-65470
28	M9-25328	M9-30328	M9-40228	M9-40328	75			M9-65275	M9-65375	M9-65475
29	M9-25329	M9-30329	M9-40229	M9-40329	80			M9-65280	M9-65380	M9-65480
30	M9-25330	M9-30330	M9-40230	M9-40330	85			M9-65285	M9-65385	M9-65485
32		M9-30332	M9-40232	M9-40332	90			M9-65290	M9-65390	M9-65490
34		M9-30334	M9-40234	M9-40334	95			M9-65295		
36		M9-30336	M9-40236	M9-40336	100			M9-652100		
38		M9-30338	M9-40238	M9-40338	105			M9-652105		
40		M9-30340	M9-40240	M9-40340	110			M9-652110		
42			M9-40242	M9-40342	115			M9-652115		
44			M9-40244	M9-40344	120			M9-652120		
45			M9-40245	M9-40345						
46			M9-40246	M9-40346						
48			M9-40248	M9-40348						
50			M9-40250	M9-40350						

Small Cannulated Screw Set

2.5mm, 3.0mm & 4.0mm

Large Cannulated Screw Set

4.5mm & 6.5mm

Thread Length

Short Thread = 25% Total Length

Long Thread = 40% Total Length

Full Thread = Total Length - 10mm

Material

Titanium Ti6Al4V



DC-0066 (CSSH-A-001)



Corporate Headquarters

1565 N. Central Expressway Suite 220

Richardson, TX 75080

P. 469.862.3030 | F. 469.862.3035 | FUSEMEDICAL.COM | OTC: FZMD