



Compression Screw Set Cannulated - Headless



- Titanium
- Self drilling
- Self tapping
- StabNGrab

MADE IN AMERICA



Compression Screw Set

Cannulated - Headless

Small Set - Implant Sizes

Length (mm)	2.0mm Headless Cannulated Screw	2.5mm Headless Cannulated Screw	3.0mm Headless Cannulated Screw
	Long Thread	Long Thread	Long Thread
10	M1-20310	M1-25310	M1-30310
11	M1-20311	M1-25311	M1-30311
12	M1-20312	M1-25312	M1-30312
13	M1-20313	M1-25313	M1-30313
14	M1-20314	M1-25314	M1-30314
15	M1-20315	M1-25315	M1-30315
16	M1-20316	M1-25316	M1-30316
17	M1-20317	M1-25317	M1-30317
18	M1-20318	M1-25318	M1-30318
19	M1-20319	M1-25319	M1-30319
20	M1-20320	M1-25320	M1-30320
21	M1-20321	M1-25321	M1-30321
22	M1-20322	M1-25322	M1-30322
23	M1-20323	M1-25323	M1-30323
24	M1-20324	M1-25324	M1-30324
25	M1-20325	M1-25325	M1-30325
26	M1-20326	M1-25326	M1-30326
27	M1-20327	M1-25327	M1-30327
28	M1-20328	M1-25328	M1-30328
29	M1-20329	M1-25329	M1-30329
30	M1-20330	M1-25330	M1-30330
32	M1-20332	M1-25332	M1-30332
34	M1-20334	M1-25334	M1-30334
36	M1-20336	M1-25336	M1-30336
38	M1-20338	M1-25338	M1-30338
40	M1-20340	M1-25340	M1-30340
45	M1-20345		
50	M1-20350		

Large Set - Implant Sizes

Length (mm)	4.5mm Headless Cannulated Screw		6.5mm Headless Cannulated Screw		
	Short Thread	Long Thread	Short Thread	Long Thread	Full Thread
30	M1-45230	M1-45330	M1-65230	M1-65330	M1-65430
32	M1-45232	M1-45332	M1-65232	M1-65332	M1-65432
34	M1-45234	M1-45334	M1-65234	M1-65334	M1-65434
36	M1-45236	M1-45336	M1-65236	M1-65336	M1-65436
38	M1-45238	M1-45338	M1-65238	M1-65338	M1-65438
40	M1-45240	M1-45340	M1-65240	M1-65340	M1-65440
42	M1-45242	M1-45342	M1-65242	M1-65342	M1-65442
44	M1-45244	M1-45344	M1-65244	M1-65344	M1-65444
46	M1-45246	M1-45346	M1-65246	M1-65346	M1-65446
48	M1-45248	M1-45348	M1-65248	M1-65348	M1-65448
50	M1-45250	M1-45350	M1-65250	M1-65350	M1-65450
52	M1-45252	M1-45352	M1-65252	M1-65352	M1-65452
54	M1-45254	M1-45354	M1-65254	M1-65354	M1-65454
56	M1-45256	M1-45356	M1-65256	M1-65356	M1-65456
58	M1-45258	M1-45358	M1-65258	M1-65358	M1-65458
60	M1-45260	M1-45360	M1-65260	M1-65360	M1-65460
65	M1-45265	M1-45365	M1-65265	M1-65365	M1-65465
70	M1-45270	M1-45370	M1-65270	M1-65370	M1-65470
75			M1-65275	M1-65375	M1-65475
80			M1-65280	M1-65380	M1-65480
85			M1-65285	M1-65385	M1-65485
90			M1-65290	M1-65390	M1-65490
95			M1-65295		
100			M1-652100		
105			M1-652105		
110			M1-652110		
115			M1-652115		
120			M1-652120		

Small Cannulated Screw Set

2.0mm & 2.5mm & 3.0mm

Large Cannulated Screw Set

4.5mm & 6.5mm

Thread Length

Short Thread = 25% Total Length

Long Thread = 40% Total Length

Full Thread = Total Length - 10mm

Material

Titanium Ti6Al4V



DC-0063 (CSSHL-A-001)



Corporate Headquarters

1565 N Central Expressway Suite 220
Richardson, TX 75080

P. 469.862.3030 | F. 469.862.3035 | FUSEMEDICAL.COM | OTC: FZMD