



medical

SUTURE ANCHOR DELIVERY SYSTEMS



GALEN


Medial Suture Anchor
Small Joint Anchor

KOPIS

Knotless Anchor

VIDA

Interference Screws



Ti6Al4V ELI
ZENIVA™ PEEK
MADE IN AMERICA



medical

SUTURE ANCHOR DELIVERY SYSTEMS

GALEN SUTURE ANCHOR

- Medial and Small Joint suture anchors in double and triple loaded platforms
- Available in PEEK and titanium with multiple needle options
- Sizes 4.5, 5.5, and 6.5 (mm) diameters
- PEEK labral push-in anchor offered in 2.9 and 3.3 (mm) diameters

KOPIS KNOTLESS ANCHOR

- Double eyelet and two-piece twist in design
- Available in both PEEK and titanium
- Sizes 4.5, 5.5, and 6.5 (mm) diameters

VIDA PEEK INTERFERENCE SCREWS

- Available in sizes 7-12 (mm) diameter and 23-30 (mm) lengths
- Mini screws available in sizes 4.5-8 (mm) diameters and 12-15 (mm) lengths

Corporate Headquarters

Fuse Medical Inc.

1565 N Central Expressway Suite 220

Richardson, TX 75080

P. 469.862.3030 | F. 469.862.3035 | FUSEMEDICAL.COM | OTC: FZMD