



ACIF CERVICAL INTERBODY SPACER


TYPEEK™ TITANIUM PLASMA SPRAY

PEEK-OPTIMA™ CORE


BIOACTIVE SURFACES*

OPTIMIZED TOOTH PROFILE

LARGE GRAFT WINDOW




MADE IN AMERICA



*Titanium plasma sprayed coating on PEEK has been shown to encourage osseointegration [Walsh et. al., 2015]


*Titanium is "Inherently Bioactive" [D.J. Hickey et al., 2015]

All TyPEEK™ implants are packaged sterile: Shaded items are special order



Height	D x (D-2.5) - Lordotic 5°	
	14 x 11.5	16 x 13.5
5	M2-C14555-S	M2-C16555-S
6	M2-C14655-S	M2-C16655-S
7	M2-C14755-S	M2-C16755-S
8	M2-C14855-S	M2-C16855-S
9	M2-C14955-S	M2-C16955-S
10	M2-C141055-S	M2-C161055-S
11	M2-C141155-S	M2-C161155-S
12	M2-C141255-S	M2-C161255-S

Height	D x (D-2.5) - Parallel (0°)	
	14 x 11.5	16 x 13.5
5	M2-C14505-S	M2-C16505-S
6	M2-C14605-S	M2-C16605-S
7	M2-C14705-S	M2-C16705-S
8	M2-C14805-S	M2-C16805-S
9	M2-C14905-S	M2-C16905-S
10	M2-C141005-S	M2-C161005-S
11	M2-C141105-S	M2-C161105-S



Height	Approximate Autograft Volume
	14 x 11.5
5	0.22 cc
6	0.27 cc
7	0.31 cc
8	0.36 cc
9	0.40 cc
10	0.45 cc
11	0.49 cc
12	0.54 cc

Corporate Headquarters
Fuse Medical Inc.

1565 N Central Expressway Suite 220
Richardson, TX 75080