



TLIF INTERBODY SPACER

TPLIF INTERBODY SPACER


TYPEEK™ TITANIUM PLASMA SPRAY*


PEEK-OPTIMA™ CORE

BIOACTIVE SURFACES**

LARGE GRAFT WINDOW

OPTIMIZED TOOTH PROFILE





Interbody System - TLIF

Height	L x W - Angle 5°
	30 x 10
7mm	M2-T307502-S
8mm	M2-T308502-S
9mm	M2-T309502-S
10mm	M2-T3010502-S
11mm	M2-T3011502-S
12mm	M2-T3012502-S
13mm	M2-T3013502-S
14mm	M2-T3014502-S
15mm	M2-T3015502-S

Interbody System - TPLIF

Height	L x W - Angle 0°	L x W - Angle 8°	L x W - Angle 0°
	28 x 10	28 x 10	32 x 10
7mm	M2-P287410-S	M2-P287810-S	M2-P327010-S
8mm	M2-P288410-S	M2-P288810-S	M2-P328010-S
9mm	M2-P289410-S	M2-P289810-S	M2-P329010-S
10mm	M2-P2810410-S	M2-P2810810-S	M2-P321001-S
11mm	M2-P2811410-S	M2-P2811810-S	M2-P3211010-S
12mm	M2-P2812410-S	M2-P2812810-S	M2-P3212010-S
13mm	M2-P2813410-S	M2-P2813810-S	M2-P3213010-S
14mm	M2-P2814410-S	M2-P2814810-S	M2-P3214010-S
15mm	M2-P2815410-S	M2-P2815810-S	M2-P3215010-S

Height	TLIF Approximate Autograft Volume (L x W)
	30 x 10mm
7mm	0.69 cc
8mm	0.80 cc
9mm	0.90 cc
10mm	1.00 cc
11mm	1.10 cc
12mm	1.21 cc
13mm	1.31 cc
14mm	1.41 cc
15mm	1.51 cc

Height	TPLIF Approximate Autograft Volume (LxW)	
	28 x 10mm	32 x 10mm
7mm	0.7 cc	0.8 cc
8mm	0.8 cc	1.0 cc
9mm	0.9 cc	1.1 cc
10mm	1.0 cc	1.2 cc
11mm	1.1 cc	1.3 cc
12mm	1.2 cc	1.5 cc
13mm	1.3 cc	1.6 cc
14mm	1.4 cc	1.7 cc
15mm	1.5 cc	1.8 cc

*Titanium plasma sprayed coating on PEEK has been shown to encourage osseointegration [Walsh et. al., 2015]

**Titanium is "Inherently Bioactive" [D.J. Hickey et al., 2015]

Corporate Headquarters
Fuse Medical Inc.

1565 N Central Expressway Suite 220
Richardson, TX 75080