


fuse ACP



Low Profile
Large Graft Window
1-5 Level Plate Options
Dual Screw Locking Mechanism




Features & Benefits

Large windows provide visualization for both the disc space and vertebral endplate.


One step, dual screw locking mechanism to prevent backout.

Low-profile, smooth, tapered edges.



Screw Angulation

4.0mm Variable Screw



4.0mm Fixed Screw




Plate Dimensions

Lengths :

1 Level : 08mm - 26mm
2 Level : 22mm - 46mm
3 Level : 36mm - 69mm
4 Level : 46mm - 78mm
5 Level : 71mm - 91mm

Radius Curvature :

Sagittal Plane
1 Level : 120mm
2 - 5 Levels : 200mm
Axial Plane : 25mm

Width :

Screw Holes: 16.5mm
Neck: 14mm

Thickness : 2.0mm


Screw Dimensions

Diameters: 4.0mm and 4.35mm

Lengths: 12mm - 18mm
(2mm increments, hole to hole)

Dual Lead Design

Color Coded Screws



Fuse Medical, Inc.

1565 N. Central Expressway, Suite 220
Richardson, TX 75080
469.862.3030

info@fusemedical.com
www.fusemedical.com