




PSS

PEDICLE SCREW SYSTEM



OPEN PEDICLE SCREW SYSTEM



PSS

PEDICLE SCREW SYSTEM

Implant Features




Tall-top percutaneous

- Polyaxial design
- Dual-lead thread
- Self-tapping
- Cannulated
- 20mm internal reduction
- Breakaway tabs

Tall-Top Screw	
DIAMETER	LENGTH
5.5	35-50mm
6.5	30-60mm
7.5	30-60mm

MIS Rods	
DIAMETER	LENGTH
5.5	40-150mm
	*5mm increments



Torx 25 driver to avoid stripping



Square thread design to prevent cross-threading



MIS Rod Inserter