




# J-PASS<sup>TM</sup> & M-FIX<sup>TM</sup>



CORACOID SOLUTIONS




Acromioclavicular Joint Repair System

## System Highlights

- Streamlined Sub-Coracoid Suture Passage
- Woven Suture Band for Stress Riser Reduction
- No Bone Tunnels Required
- Improved Biomechanics



The **J-PASS™** unique suture delivery system provides reproducible **270 degree** suture-passing around the coracoid and clavicle.



The **M-FIX™** provides broad, stable fixation with **NO DRILL HOLES** and offers compatibility with graft and fracture plates.

Provisional fixation to facilitate adjustments prior to final fixation.

FG1001	J-Pass™ Passing Instrument
FG3004	4mm M-Fix™ Implant & Accessories
FG3007	7mm M-Fix™ Implant & Accessories

**Fuse Medical, Inc.**

1565 N. Central Expressway, Suite 220  
Richardson, TX 75080  
469.862.3030

info@fusemedical.com  
www.fusemedical.com

JPMF-B-001 REV A Effective date 7/14/2022 DC-0230