



POSTERIOR LUMBAR INTERBODY FUSION


TYPEEK™ TITANIUM PLASMA SPRAY

ANATOMICAL ENDPLATE CONTACTS


SIMPLE INSTRUMENTATION

PEEK-OPTIMA™ CORE

LARGE GRAFT WINDOW



Titanium Plasma Spray*



Interbody System - PLIF


Height	L x W - Angle 0°	L x W - Angle 8°
	22 x 10mm	22 x 10mm
7mm	M02-P227010-S	M02-P227810-S
8mm	M02-P228010-S	M02-P228810-S
9mm	M02-P229010-S	M02-P229810-S
10mm	M02-P2210010-S	M02-P2210810-S
11mm	M02-P2211010-S	M02-P2211810-S
12mm	M02-P2212010-S	M02-P2212810-S
13mm	M02-P2213010-S	M02-P2213810-S
14mm	M02-P2214010-S	M02-P2214810-S
15mm	M02-P2215010-S	M02-P2215810-S

Interbody System - PLIF

Height	L x W - Angle 0°	L x W - Angle 4°	L x W - Angle 8
	24 x 10mm	24 x 10mm	24 x 10mm
7mm	M02-P247010-S	M02-P247410-S	M02-P247810-S
8mm	M02-P248010-S	M02-P248410-S	M02-P248810-S
9mm	M02-P249010-S	M02-P249410-S	M02-P249810-S
10mm	M02-P2410010-S	M02-P2410410-S	M02-P2410810-S
11mm	M02-P2411010-S	M02-P2411410-S	M02-P2411810-S
12mm	M02-P2412010-S	M02-P2412410-S	M02-P2412810-S
13mm	M02-P2413010-S	M02-P2413410-S	M02-P2413810-S
14mm	M02-P2414010-S	M02-P2414410-S	M02-P2414810-S
15mm	M02-P2415010-S	M02-P2415410-S	M02-P2415810-S

Interbody System - PLIF

Height	L x W - Angle 0°	L x W - Angle 8°
	26 x 10mm	26 x 10mm
7mm	M02-P267010-S	M02-P267810-S
8mm	M02-P268010-S	M02-P268810-S
9mm	M02-P269010-S	M02-P269810-S
10mm	M02-P2610010-S	M02-P2610810-S
11mm	M02-P2611010-S	M02-P2611810-S
12mm	M02-P2612010-S	M02-P2612810-S
13mm	M02-P2613010-S	M02-P2613810-S
14mm	M02-P2614010-S	M02-P2614810-S
15mm	M02-P2615010-S	M02-P2615810-S



*Titanium plasma sprayed coating on PEEK has been shown to encourage osseointegration [Walsh et. al., 2015]

**Titanium is "Inherently Bioactive" [D.J. Hickey et al., 2015]