



DEMINERALIZED BONE MATRIX

FusePure™ Demineralized Bone Matrix Biologic scaffold is composed of 100% pure human bone.

- All-natural demineralized tissue with no synthetic carriers to direct and support bone formation
- Patient safety ensured through terminal sterilization (SAL 10^{-6})
- Stored at room temperature for convenient storage and distribution
- Two-year shelf life



FusePure - DBM Putty

- Cohesive & conformable for reliable implantation
- Easily mixes with fluid or allograft/autograft for a variety of applications
- Ready to use

FusePure Plus - DBM Putty with Cancellous Chips

- Added cancellous bone provides an osteoconductive scaffold
- Easily mixes with fluid or allograft/autograft for a variety of applications
- Ready to use

FusePure Fiber - DBM Putty with Cortical Fibers

- Fibrous connection provides superior handling over traditional putty
- Alleviates washout and/or migration
- Ready to use
- Easily mixes with fluid or allograft/autograft for a variety of applications



FusePure™ - DBM Putty

	PRODUCT NUMBER	SIZE
FusePure™ - DBM Putty	FZPD-001	1.0cc
FusePure™ - DBM Putty	FZPD-002	2.5cc
FusePure™ - DBM Putty	FZPD-005	5.0cc
FusePure™ - DBM Putty	FZPD-010	10cc

FusePure™ Plus - DBM Putty with Cancellous Chips

	PRODUCT NUMBER	SIZE
FusePure™ - DBM Putty with Cancellous	FZPP-001	1.0cc
FusePure™ - DBM Putty with Cancellous	FZPP-002	2.5cc
FusePure™ - DBM Putty with Cancellous	FZPP-005	5.0cc
FusePure™ - DBM Putty with Cancellous	FZPP-010	10cc

FusePure™ Fiber - DBM Putty with Fibers

	PRODUCT NUMBER	SIZE
FusePure™ Fiber - DBM Putty with Fibers	FZPF-001	1.0cc
FusePure™ Fiber - DBM Putty with Fibers	FZPF-002	2.5cc
FusePure™ Fiber - DBM Putty with Fibers	FZPF-005	5.0cc
FusePure™ Fiber - DBM Putty with Fibers	FZPF-010	10cc

Fuse Medical, Inc.

4343 Sigma Rd., Suite 500
Farmers Branch, TX 75244
469.862.3030

info@fusemedical.com
www.fusemedical.com

FPP-001 REV B DC-0468