




FUSEFIX HAMMER TOE

TITANIUM


STERILE PACKAGED

STRAIGHT & ANGLED OPTIONS

SOLID & CANNULATED OPTIONS



Titanium Alloy
Ti6Al4V
MADE IN AMERICA



- Built in AO Connection / disposable handle
- No additional instrumentation for a faster procedure
- Integrated instruments including drill & broach
- Barbed distal end for enhanced bone purchase & fixation
- Screws contain self-drilling and tapping features
- Breakaway design

Implant - Ø2.0 Screw

Ø2.0 Screw		
Part #	Length	Description
M 13-C20013	13mm	FuseFix Cannulated Ø2.0mm x 13mm
M 13-S20013	13mm	FuseFix Solid 0° Ø2.0mm x 13mm
M 13-S20113	13mm	FuseFix Solid 10° Ø2.0mm x 13mm

Implant - Ø2.5 Screw

Ø2.5 Screw		
Part #	Length	Description
M 13-C25015	15mm	FuseFix Cannulated Ø2.5mm x 15mm
M 13-S25015	15mm	FuseFix Solid 0° Ø2.5mm x 15mm
M 13-S25115	15mm	FuseFix Solid 10° Ø2.5mm x 15mm

Implant - Ø3.0 Screw

Ø3.0 Screw		
Part #	Length	Description
M13-C30015	15mm	FuseFix Cannulated Ø3.0mm x 15mm
M 13-S30015	15mm	FuseFix Solid 0° Ø3.0mm x 15mm
M 13-S30115	15mm	FuseFix Solid 10° Ø3.0mm x 15mm

Instruments - Sterile / Disposable	
Part #	Description
M19-01000-S	Mini Fixed AO Handle
M19-01009-S	Mini Fixed AO Handle w/ 0.9mm guide wire (6") kit
M19-01011-S	Mini Fixed AO Handle w/ 1.1mm guide wire (6") kit

Corporate Headquarters

Fuse Medical Inc.

4343 Sigma Rd., Suite 500

Farmers Branch, TX 75244