



TyWEDGE™ SYSTEM

EVANS AND COTTON PROCEDURES

TYPEEK™ TITANIUM PLASMA SPRAY





PEEK RADIOLUCENT BODY

BIOACTIVE SURFACES


OPTIMIZED TOOTH PROFILE

BONE GRAFT CONTAINMENT

COMPLETE INSTRUMENTATION




MADE IN AMERICA



*Titanium plasma sprayed coating on PEEK has been shown to encourage osseointegration [Walsh et. al., 2015]


Implants - Evans



Height	Length x Width		
	18 x 18mm	20 x 20mm	22 x 22mm
8	X06-E01808	X06-E02008	X06-E02208
10	X06-E01810	X06-E02010	X06-E02210
12	X06-E01812	X06-E02012	X06-E02212

Height	Approximate Autograft Volume		
	18 x 18mm	20 x 20mm	22 x 22mm
8	0.51 cc	0.74 cc	1.00 cc
10	0.59 cc	0.85 cc	1.15 cc
12	0.67 cc	0.96 cc	1.31 cc

Implants - Cotton



Height	Length x Width	
	16 x 14mm	20 x 14mm
4.5	X06-C01604	X06-C02004
5.5	X06-C01605	X06-C02005
6.5	X06-C01606	X06-C02006

Height	Approximate Autograft Volume	
	16 x 14mm	20 x 14mm
4.5	0.17 cc	0.25 cc
5.5	0.20 cc	0.27 cc
6.5	0.22 cc	0.31 cc

Corporate Headquarters

Fuse Medical Inc.

4343 Sigma Rd., Suite 500

Farmers Branch, TX 75244

TWS-B-001 REV C DC-0470

P. 469.862.3030 | F. 469.862.3035 | FUSEMEDICAL.COM | OTC: FZMD