




PSS

PEDICLE SCREW SYSTEM





PSS

PEDICLE SCREW SYSTEM

Implant Features



Tall-Top Percutaneous

- Polyaxial design
- Dual-lead thread
- Self-tapping
- Cannulated
- 20mm internal reduction
- Reinforced breakaway tabs

Tall-Top Screw	
DIAMETER	LENGTH
5.5	35-50mm
6.5	30-60mm
7.5	30-60mm

MIS Pre-Curved Rods	
DIAMETER	LENGTH
5.5	40-150mm
	*5mm increments



Torx 25 driver to avoid stripping

Square thread design to prevent cross-threading



MIS Rod Inserter



MIS Rod Reducer

