

Clench™ Compression Staple



Sterile - Single Use Nitinol Compression Staple System

07 X 07




09 X 07




11 X 10




15 X 12



18 X 18




20 X 20



Technical Features

- No heat activation required
- Constant compression
- Smooth profile
- All-in-one instruments & implant kit
- Sterile, single use kits
- Reloadable instrument design
- Multiple tamp options
- Robust instrumentation



SIZE	WIDTH (W) mm	LENGTH (L) mm	THICKNESS mm	COMPRESSION mm
7mm x 7mm	7mm	7mm	1.5	1.9
9mm x 7mm	9mm	7mm		1.9
11mm x 10mm	11mm	10mm		2.8
15mm x 12mm	15mm	12mm	2.3	3.4
18mm x 18mm	18mm	18mm		5.3
20mm x 20mm	20mm	20mm		5.9







PART #	CONTENTS
152-40001	(1) Implant Sizer

PART #	SIZE	CONTENTS
152-30707	7mm x 7mm	(1) Drill Guide (1) Drill
152-30907	9mm x 7mm	(2) Guide Post (1) Delivery Instrument
152-31110	11mm x 10mm	(1) Nitinol Staple
152-31512	15mm x 12mm	Note: Contents are specific to the implant size included. Do not use instruments with other implant sizes.
152-31818	18mm x 18mm	
152-32020	20mm x 20mm	

PART #	SIZE	CONTENTS
152-40707	7mm x 7mm	(1) Delivery Instrument (1) Nitinol Staple
152-40907	9mm x 7mm	Note: Contents are specific to the implant size included. Do not use instruments with other implant sizes.
152-41110	11mm x 10mm	
152-41512	15mm x 12mm	
152-41818	18mm x 18mm	
152-42020	20mm x 20mm	

SURGICAL PROCEDURE



This brochure contains information intended to highlight a Reign Medical product. It is highly recommended to review the package insert, surgical technique guide, and/or participate in a product training course before using any Reign Medical product. Surgeons should use their medical experience when choosing which products and techniques to use for each patient. Please contact your authorized representative for any questions and product availability.



Distributed By:

Fuse Medical Inc
4343 Sigma Road, Suite 500
Farmers Branch, TX 75244
BRO-CCS-001-REV A DC-0513

www.fusemedical.com
(O) 469-862-3030
(F) 469-862-3035

Manufactured By:

Reign Medical www.reignmedical.com
6303 E. 102nd Street (O) 918-302-9546
Tulsa, OK 74137 (F) 918-844-4356