



COMPRESSION SCREW SYSTEMS

HEADED AND HEADLESS CANNULATED SCREWS



- Titanium
- Self-drilling
- Self-tapping
- Pick-n-Stick



COMPRESSION SCREW SYSTEMS

HEADED AND HEADLESS CANNULATED SCREWS

Headless, Small Set		
Size	Length	Thread
2.0mm	10-50mm	Long
2.5mm	10-40mm	Long
3.0mm	10-40mm	Long

Headed, Small Set		
Size	Length	Thread
2.5mm	10-30mm	Long
3.0mm	10-40mm	Long
4.0mm	14-50mm	Short
4.0mm	14-50mm	Long

Headless, Large Set		
Size	Length	Thread
4.5mm	30-70mm	Short
4.5mm	30-70mm	Long
6.5mm	30-120mm	Short
6.5mm	30-90mm	Long
6.5mm	30-90mm	Full

Headed, Large Set		
Size	Length	Thread
4.5mm	30-70mm	Short
4.5mm	30-70mm	Long
6.5mm	30-120mm	Short
6.5mm	30-90mm	Long
6.5mm	30-90mm	Full

Thread Length

Short Thread = 25% Total Length

Long Thread = 40% Total Length

Full Thread = Total Length - 10mm

Material

Titanium Ti6Al4V

Drivers

Torx

Fuse Medical, Inc.

4343 Sigma Road, Suite 500

Farmers Branch, TX 75244

469.862.3030

BRO-CSS001

Rev A DC-0514