




# PSS

PEDICLE SCREW SYSTEM



MIS PEDICLE SCREW SYSTEM



# PSS

PEDICLE SCREW SYSTEM

## Implant Features




### Tall-Top Percutaneous

- Polyaxial design
- Dual-lead thread
- Self-tapping
- Cannulated
- 20mm internal reduction
- Reinforced breakaway tabs


Tall-Top Screw	
DIAMETER	LENGTH
5.5	35-50mm
6.5	30-60mm
7.5	30-60mm

MIS Pre-Curved Rods	
DIAMETER	LENGTH
5.5	40-150mm
	*5mm increments



Torx 25 driver to avoid stripping

Square thread design to prevent cross-threading



### MIS Rod Inserter



### MIS NAV Driver



#### Fuse Medical, Inc.

4343 Sigma Rd., Suite 500  
Farmers Branch, TX 75244  
469.862.3030

info@fusemedical.com  
www.fusemedical.com