




# PSS

PEDICLE SCREW SYSTEM




OPEN PEDICLE SCREW SYSTEM



# PSS

PEDICLE SCREW SYSTEM

## Implant Features



### Low-top Option

- Low profile
- Polyaxial design
- Dual-lead thread
- Self-tapping



### Solid and Cannulated Options

Low-Top Screw	
DIAMETER	LENGTH
4.5	25-50mm
5.5	25-55mm
6.5	30-60mm
7.5	30-60mm


Mid-Top Screw	
DIAMETER	LENGTH
4.5	25-50mm
5.5	25-55mm
6.5	30-60mm
7.5	30-60mm

Pre-Curved Rods	
DIAMETER	LENGTH
5.5	35-150mm
	*5mm increments

Straight Rods	
DIAMETER	LENGTH
5.5	150, 200-500mm
	*100mm increments

### Reduction Option

- 20mm internal
- Breakaway tabs



**Fuse Medical, Inc.**  
 4343 Sigma Rd., Suite 500  
 Farmers Branch, TX 75244  
 469.862.3030

info@fusemedical.com  
 www.fusemedical.com